MEDICAL CONTROL POLICY STATEMENT/ADVISORY



No. 2023-06 **Date:** August 8th, 2023

Naloxone Syringes

All EMS Providers:

As a result of continuing medication shortages, some facilities have been utilizing different naloxone syringes. The syringe can be used both with our atomizers and our IV tubing.

When utilizing these syringes for our nasal atomizers, the atomizer can be directly connected to the syringe without any further steps.



If you need to connect the syringe to one of the hubs of your IV tubing, or directly to the clave adaptor on your IV, you will need to remove the end of the syringe.



If you have any questions, do not hesitate to reach out to us.